



Model 175 Cot Fastener, Plus Mounting Set 175-1, 175-2, 175-3, 175-4, and 175-5

Model 175 consists of the Antler and Rail only. Mounting Sets 175-1, 175-2, 175-3, 175-4, and 175-5 consists of blocks, bracket, hardware, and floor plates necessary to fasten the antler and rail to the floor or wall of the vehicle. The Fastener Kit Complete includes all: antler, rail, blocks, hardware, and floor plates.


The Complete Kit (antler, rail, and mounting set) offers a full line of fasteners to secure a cot in an ambulance. The following features are included on all models:

- Wall and Floor mounting
- Permanent and detachable mounting
- Detachable models allow easy floor cleaning
- Wall mount models provide maximum aisle space
- Antler and Rail made of aluminum (stainless steel optional)
- Compatible with cot series Model 25, 26, 30, 35, 93, and POWERFlexx

Model 175 plus Mounting Set 175-1

- Model 175 plus Mounting Set 175-1
- Antler: Permanent Floor Mount
- Rail: Permanent Wall Mount

0561750 Model 175, Antler and Rail only
0561751 Model 175-1, Mounting Set only
0001751 175-1 Fastener Kit Complete



Model 175 plus Mounting Set 175-2

- Model 175 plus Mounting Set 175-2
- Antler: Permanent Floor Mount
- Rail: Permanent Floor Mount


0561750 Model 175, Antler and Rail only
0561752 Model 175-2, Mounting Set only
0001752 175-2 Fastener Kit Complete



Model 175 plus Mounting Set 175-3

- Model 175 plus Mounting Set 175-3
- Antler: Removable Floor Mount
- Rail: Removable Floor Mount

0561750 Model 175, Antler and Rail only
0561753 Model 175-3, Mounting Set only
0001753 175-3 Fastener Kit Complete



Model 175 plus Mounting Set 175-4

- Model 175 plus Mounting Set 175-4
- Antler: Removable Floor Mount
- Rail: Removable Floor Mount
- Includes an extra set of mounting plates to allow:
 - a) two different mounting locations (side or center)
 - b) setups for two different cot styles (example: Model 30 and Model 35)


0561750 Model 175, Antler and Rail only
0561754 Model 175-4, Mounting Set only
0001754 175-4 Fastener Kit Complete



Model 175 plus Mounting Set 175-5

- Model 175 plus Mounting Set 175-5
- Antler: Removable Floor Mount
- Rail: Permanent Wall Mount

0561750 Model 175, Antler and Rail only
0561755 Model 175-5, Mounting Set only
0001755 175-5 Fastener Kit Complete




Other Cot Fasteners



EFX-Lock

- Permanent floor mount cot fastener
- Compatible with EFX Model 5126 and 4050 Series Transporters


0573073 EFX Lock
0573074 EFX Lock MK II



Model 170 Cot Fastener

- Permanent floor mount for rear lock housing
- Permanent floor mount for antler
- Compatible with EFX and 4050 Series Transporters


0573059 Model 170



Model 5-14 Cot Fastener

- Permanent floor mount for rear fastener handle
- Permanent floor mount for front antler
- Compatible with Model 25 Cot


0573057 Model 5-14



Model 173 Cot Fastener Plus Mounting Kit 173-2

- Permanent floor mount for rear fastener hook
- Permanent floor mount for front antler
- Compatible with Model 28 Chair Cot

0561730 Model 173
0561731 Model 173-2, Mounting Set only



Model 173 Cot Fastener Plus Mounting Kit 173-3

- Detachable floor mount for rear fastener hook
- Detachable floor mount for front antler
- Compatible with Model 28 Chair Cot

0561730 Model 173
0561732 Model 173-3, Mounting Set only

