



## Ferno Oxy-Clip2™ Spring Loaded Oxygen Bottle Holder

- Works with the EZ Glide™ Stair Chair
- Made of lightweight aluminum with a nylon webbing and zinc die cast/nickel plated D-ring fasteners
- Spring-loaded mechanism for quick attachment and removal
- Will hold a D-size oxygen bottle in place during transport freeing up the operator's hands improving both safety and patient care

**0821977** Ferno OxyClip2™ for 59-T and 59 EZ Glide



## Ferno X-O<sub>2</sub> Oxygen Mounting System

- Works with a 35-P
- The X-O<sub>2</sub> System can utilize "D"- or "E"-size oxygen bottles. The patient surface-mount control center allows for either a low flow (with flow control from 0 to 25 LPM) or high flow connection in support of patient care

**CA5320** Ferno X-O<sub>2</sub> for 35-P



*Surface mount flow control*

## Ferno Oxygen Cylinder Holder

- Works with either a PROFlexx 35-P or 93-P
- Can carry a "D" size, "Jumbo-D", or a FernoO<sub>2</sub> DD-lite
- FernoO<sub>2</sub> DD-lite is also compatible

**0819958** for 35-P

**0821947** for 93-P



## Ferno Oxy-Clip™ Cylinder Holder

- Fits on the cot side rail (19mm diameter)
- Securely holds either steel or aluminum "D" size cylinders
- Spring-loaded latch enables the Oxy-Clip™ to instantly attach to the cot rail
- When transferring the patient to the hospital bed, the Oxy-Clip™ simply rotates down and out of the way as the cot side rail is lowered
- Will not accept Jumbo "D" size cylinder

**0085500** Model 514 for models 25, 26, 35, 93, 4140, 4150, EFX-1, 2, and 3

**0085520** Model 514-1 for model 28



## Ferno Model 274 PacRac™



- Quick-clip for fast and secure attachment to cot's frame
- Multiple mounting options from mid-section to foot-end of cot
- Nylon security straps with hook-and-loop closures keep equipment secure during patient transport
- Lightweight aluminum construction
- Weight capacity 75 lb / 34 kg

**0818933** Model 274 for 93-P, 35-P, and POWERFlexx

**0818952** Model 274-S for 28



## Ferno Model 513 Series IV Poles

- Designed for bagged or bottled IV solutions
- Elastic Velcro® strap limits sway of solution containers during transport
- Height control tab adjusts for appropriate solution levels
- Folds along cot mattress for easy storage

**Model 513-EFX** is Compatible with EFX-1, EFX-2, and EFX-3 Stretchers

**0832024** EFX

**Model 513** is compatible with Models 25, 26, 30, 35, 93 & POWERFlexx Cot:

**0085130** Model 513

**Model 513-4** is specifically made for the Model 71 Series Basket Stretcher

**0085480** Model 513-4 for Model 71 Basket Stretcher



**Model 513**



**Model 513-4**

## Ferno Model 513-10 Series Heavy Duty IV Poles

- Permanently mounts to cot's main frame
- Adjustable height for placement of IV solutions at appropriate levels
- Poles are nonmagnetic for MRI use (excluding Model 513-LBS)
- Capable of holding four one-liter containers of fluid and two transport IV pumps
- Large hook-and-loop straps stabilize IV containers
- Folds along cot mattress for easy storage (excluding Model 513-13)
- Lock sleeve secures pole in upright position

**Model 513-10** is compatible with Models 25, 26, 30, 35, 93 and POWERFlexx Cots. Length 33"-49".

**0087156** Model 513-10, HD NM

**Model 513-11** is specifically made for Ferno-Flex™ Model 28 Cot. Mounts to rectangular main frame.

**0087157** Ferno-Flex™ Model 513-11, HD NM

**Model 513-13** is compatible with Models 25, 26, 30, 35, 35-P, 93, 93-P and POWERFlexx Cots. Length 24"-38".

**0087172** Model 513-13

**Model 513-LBS** is designed for use with the LBS.

**0087173** Model 513-LBS



**Model 513-10**

**Model 513-11**

**Model 513-13**



## Ferno Model 417 Shoulder Torso Harness Restraint

- Four individually adjustable sections
- Quick-release chrome-plated metal buckle
- Should be used in conjunction with two Ferno Model 430 Patient Restraints
- Available in burgundy, orange and black

### Model 417-1

**0313915** Black

**0313896** Burgundy

**0313914** Orange

**0313996** For EFX Stretchers, Black



## Ferno Model 430 Series Cot Restraint Straps Two Piece

- Two-piece primary-use patient restraints made of 2" nylon webbing
- Quick-release automotive-type
- Chrome-plated male and black painted female metal buckle
- Resists rust and corrosion
- Push button release
- Hemmed tab on female end
- Adjustable on male section
- Available in black, burgundy, and orange

	<b>Black</b>	<b>Burgundy</b>	<b>Orange</b>
<b>5'</b>	0313892	0313890	0313891
<b>7'</b>	0313911	0313913	0313912
<b>9'</b>	0313305	0313306	0312560



## Ferno 414-OL Over-The-Lap Wrist Restraints

- Adjustable to hold patient's wrists securely in lap
- Restraint fastens to both sides of the cot
- Hook-and-loop closures
- Available in orange and burgundy

**0312320** Model 414-OL Orange

**0313090** Model 414-OL Burgandy



## Ferno 414-R Wrist Restraints, Cot Sidearm

- Restraints hold the patient's wrists to the top of the cot's side rails
- Hook-and-loop closures
- Set of two available in orange and burgundy

**0312310** Model 414-R, set of two, Orange

**0313088** Model 414-R, set of two, Burgundy



## Ferno 415 Ankle Restraints

- Fastens to both sides of the cot's main frame
- Holds patient's ankles securely in place
- Restraints allow each leg to be individually secured or released

**0312330** Model 415-AR Orange

**0313089** Model 415-AR Burgundy





# Ferno Blankets & Pillow

## Model 352 100% Acrylic-Blend Blanket

- Fiberwoven, material (acrylic) is stitched around perimeter to withstand repeated washings
  - Measures 72" X 90"
  - Machine wash and dry
- 0311140** Model 352, Blue only

352



## Model 353 Quilted Nylon Blanket

- 100% nylon foam and 100% cotton broadcloth
  - Machine wash and dry or dry clean
  - Measures 40" x 90"
- 0311150** Model 353, Red only

353



## Model 354 Wool-Blend Blanket

- Lightweight: 70% wool, 30% synthetic fibers
  - Dry clean preferred or machine wash, air dry
  - Fire retardant to CS 191-53
  - Measures 60" X 80"
- 0311180** Model 354, Charcoal Gray only

354



## Model 355 100% Wool Blanket

- High-quality wool with surged perimeter provides excellent insulation.
  - Dry clean preferred or machine wash, hang dry
  - Meets Class 1 Fire Rating naturally
  - Measures 50" x 90"
- 0311190** Model 355, Red only

355



## Model 356 Fire Blanket Only

- 70% wool, 30% synthetic fibers
  - Dry clean preferred or machine wash, hang dry
  - Fire retardant (CS191-53)
  - Also available with wall case
  - Measures 60" X 84"
- 0311200** Model 356, Fire Blanket, Burgundy

## Model 680 Fire Blanket with Wall Case

**0312740** Model 680, Fire Blanket (burg) with Orange Case



Model 680  
blanket and wall case

Model 356  
blanket only

## Model 366 Cot Pillow

- Durable, waterproof cover encloses soft, synthetic filler
- Fits into a pillowcase sideways, leaving a generous flap that tucks under the mattress and secures the pillow
- Nonallergenic

366

**0311310** Model 366, White only

