



Ferno Model 59-T EZ-Glide™ Evacuation Stair Chair (with tracks)

The EZ Glide chair can be used as a stair chair, evacuation chair, and transport chair.

- Two rear tracks contact the lip of stairs allowing them to glide over the edge of each step eliminating the need to lift chair and patient
- All metal frame construction allows you to carry bariatric patients
- Available in five colors (blue, red, orange, green, and platinum)
- Two 4" front swivel wheels and two 6" rear wheels, with locks
- Five position rear lift bar (from 37.5" to 63.5" / 95 - 161 cm)
- Five position front telescoping handles (from 5.4" to 19" / 14 - 48 cm)
- Seat area available in two versions, snap-on ABS plastic panels (shown) or a removable one-piece vinyl cover
- Foldable footrest



<i>EZ-Glide with Tracks</i>	ABS Panel	Vinyl Cover
<i>59-T Blue with restraints</i>	00059TB	
<i>59-T Red</i>	0731305	0731310
<i>59-T Orange</i>	0731306	0731311
<i>59-T Green</i>	0731307	0731312
<i>59-T Blue</i>	0731308	0731313
<i>59-T Red with restraints</i>	00059TR	
<i>59-T Platinum</i>	0731309	0731314


	Height	Width	Max Depth	Min Depth	Seat Depth	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping L,W,H
<i>Imperial</i>	37½	20.3	51	28	16	8	33	500	41	40x21x9
<i>Metric</i>	95	52	130	73	41	20	15	227	19	102x53x23

Kwik-Klip Restr.



P/N 0314002

430P Restraints



0313797

IV Bag Holder



0821976

Oxygen Holder



0821977

Headrest



0314000

Non-Lock Handle



0821974

Lock Handle



0821975

Model 59 EZ-Glide™ Evacuation Stair Chair (without tracks)



The EZ Glide chair can be used as a stair chair, evacuation chair, and transport chair.

- Two rear tracks contact the lip of stairs allowing them to glide over the edge of each step eliminating the need to lift chair and patient
- All metal frame construction allows you to carry bariatric patients
- Available in five colors (blue, red, orange, green, and platinum)
- Two 4" front swivel wheels and two 6" rear wheels, with locks
- Five position rear lift bar (37.5" to 63.5" / 95 - 161 cm) for operators of various heights
- Five position front telescoping handles (from 5.4" to 19" / 14 - 48 cm)
- Seat area available in two versions, snap-on ABS plastic panels (shown) or a removable one-piece vinyl cover
- Foldable footrest



0822072 Wall Mount Bracket



<i>EZ-Glide without Tracks</i>	ABS Panel	Vinyl Cover
59 Red	0731315	0731320
59 Orange	0731316	0731321
59 Green	0731317	0731322
59 Blue	0731318	0731323
59 Platinum	0731319	0731324


	Height	Width	Max Depth	Min Depth	Seat Dept	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	37½	20.3	40	27	16	8	26	500	41	40x21x9
Metric	95	52	102	69	41	20	12	227	19	102x53x23

Kwik-Klip Restr.



P/N 0314002

430P Restraints



0313797

IV Bag Holder



0821976

Oxygen Holder



0821977

Headrest



0314000

Non-Lock Handle



0821974

Lock Handle



0821975



Ferno Model 40 Stair Chairs

- Easy to maneuver by two EMTs in-and-out of cramped quarters
- Two 5" stationary rear wheels
- Telescoping front handles with two-position
- Vinyl-coated nylon cover is simple to clean, resists stains, blood and bodily fluids
- Three 7' quick-release patient restraints hold the patient securely
- Chair mounts to interior or exterior walls
- Ideal for high-rise buildings or industrial settings
- Available molded ABS plastic seat and back or with the original soft vinyl cover
- Available in burgundy and orange

0714000 Model 40 with Vinyl Cover, Orange

0714856 Model 40-WP with ABS Panels, Orange

0714801 Model 40 with Vinyl Cover, Burgundy



Model 40 with durable
ABS Panels



Model 40 with Vinyl Cover

Telescoping handles enable
front-end EMT to move
without interfering with the
patient's feet



	Height	Width	Open Depth	Folded Depth	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	36	20	27	7	19	350	24	38x21x8
Metric	91	51	69	18	9	159	11	97x53x20

Model 40-OS

- Economy version of the Model 40
- Without telescoping front handles
- Two 5" stationary rear wheels
- Four lifting handles
- Locking safety latch
- Vinyl covers only

0714841 Model 40-OS with Vinyl Cover, Orange



Model 40-OS



	Width	Height	Open Depth	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	20	36	27	7	14	350	21	38x21x8
Metric	51	91	69	18	9	159	10	97x21x8

Ferno Model 42/48 Stair Chairs



- Designed to be stored in small places
- Functions well in confined areas such as narrow hallways and multiple-landing stairways
- Two 4" rear stationary wheels and two 2" low-profile front guide wheels
- Roll over most carpeted or hard surfaces
- Four lifting handles with nonslip handgrips
- Soft vinyl-coated cover is simple and easy to clean, resists blood, stains and bodily fluids
- Features a locking safety latch
- Safety latch is manually released for storage
- Two 7' quick-release patient restraints securely hold the patient

Model 42 is available with a durable, easy-to-clean, molded ABS plastic seat and back, with built-in clipboard holder, or with the original soft vinyl cover. Available in burgundy and orange.

- 0714200** Model 42 with Vinyl Cover, Orange
- 0714802** Model 42 with Vinyl Cover, Burgundy
- 0714857** Model 42-WP with ABS Panels, Orange
- 0714859** Model 42-WP with ABS Panels, Burgundy



	Height	Width	Open Depth	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	37	19	24	6	19	350	23	26x20x7
Metric	95	49	61	16	9	159	10	66x51x18

Model 48 Stair Chair is identical to Model 42, only 3" narrower and available with soft vinyl cover only. Ideal for use in airplanes, boats, restaurants and all areas where access is severely limited. Available in burgundy and orange.

- 0714800** Model 48 with Vinyl Cover, Orange
- 0714803** Model 48 with Vinyl Cover, Burgundy



	Height	Width	Open Depth	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	37	16	24	6	18	350	23	26x20x7
Metric	95	40¼	61	16	8	159	10	66x51x18



Ferno Model 49 Sirocco™ Evacuation Chair

- Innovative tri-wheel assembly helps to reduce operator strain
- Can be used to ascend and descend stairs by rolling instead of carrying a patient
- Wheel retarding system enables the operator to control the speed while descending stairs
- Foam grip handles on foot-end and head-end
- Strapping system, adjustable shoulder harness, lap strap and integral foot restraint are included as standard
- Optional ratcheting headrest and 3-point harness restraint system available

0714861 Model 49 Sirocco Chair without headrest

0819955 Optional Ratcheting Headrest

0313921 Optional 3-Point Harness Restraint


0313986 Optional Chair Cover, Burgundy

Chair Cover



P/N 0313986

Ratcheting Headrest



0819955



Shown with optional ratcheting headrest and standard restraints



	Height	Width	Open Depth	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	54	20	37	15	40	350	46	45x21x24
Metric	137	51	94	38	18	159	21	114x53x61

Compact III Chair

The Compact III Chair has been designed to transport a patient, in a seated position, up or down stairs and also on level ground.

- Innovative tri-wheel system enables smooth operation up and down stairs and reduces requirement for lifting
- Upper carrying handle is fully adjustable allowing easy operation and helps improve operator control

0714853 Compact III Chair



	Height	Stored Height	Width	Open Depth	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping H,L,W
Imperial	39	36	16	17	43	28	280	46	45x2x24
Metric	98	92	41	43	36	13	127	21	45x2x24