

Ferno Model 12 Emergency Stretcher



- Exclusive slimline center hinges enable stretcher to fold for compact storage
- Heavy duty vinyl-coated nylon cover
- Lightweight aluminum frame
- Can be used with Model 385-A Storage/Carrying Case
- Includes two 5' quick-release patient restraints



CE

0107987 Model 12, Burgundy **0311650** Storage Case, Burgundy
0107986 Model 12, Orange **0311630** Storage Case, Orange

	Open	Folded	Width	Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	74	37	19	1	9	350	12	38x21x4
Metric	187	94	47	3	4	159	5	97x53x1

Ferno Model 12-C Emergency Breakaway Stretcher

- Low profile frame for easier log-rolling or lifting a patient
- Stretcher halves separate for easy removal without moving the patient
- Frame adjusts to align with ambulance cot backrest enabling patient to have head elevated during transport
- Includes two 5' quick-release patient restraints
- Heavy duty vinyl-coated nylon cover



CE

0107743 Model 12-C, Burgundy
0107744 Model 12-C, Orange

	Open Length	Folded Length	Width	Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	74	37	19	9	15	350	15	38x21x4
Metric	187	94	47	4	7	159	7	97x53x10

Ferno Model 11 Emergency Stretcher with Wheels

- Convenient trigger-release
- Folding 4" wheels and posts can be used on the street and then placed directly onto a squad bench equipped with Model 533 Wheel Cups and Models 534/535 Post Cups
- Includes two 5' quick-release patient restraints
- Heavy duty vinyl-coated nylon cover
- Can be used with Model 385-A Storage/Carrying Case




CE

0107985 Model 11, Burgundy **0311650** Storage Case, Burgundy
0107984 Model 11, Orange **0311630** Storage Case, Orange

	Open Length	Folded Length	Width	Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	74	37	19	8	15	350	18	38x21x4
Metric	187	94	47	21	7	159	8	97x53x10

Ferno Model 9 Adjustable Backrest Emergency Stretcher with Wheels

- Head-end frame section telescopes to shorten stretcher
- Backrest is adjustable to eight angles (0° to 75°)
- Four-inch wheels at head-end
- Stationary posts at foot-end
- Folds in half for compact storage
- Can be used with Model 385-A Storage/Carrying Case
- Includes two 5' quick-release patient restraints
- Heavy duty vinyl-coated nylon cover



Adjustable Backrest

CE

0107731 Model 9, Burgundy
0100900 Model 9, Orange

	Open Length	Folded Length	Width	Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	74	371	19	6	17	350	18	41x20x8
Metric	187	94	47	16	8	159	8	104x51x20



Ferno 107 Series Combination Stretcher Chairs

- Vinyl-coated nylon cover
- Folds in half for compact storage
- High-back chair design
- Folding handles at foot-end provide excellent grip
- Stationary or hinged 4" wheels
- Folding head-end posts double as lift handles and support legs
- Legs fold underneath the stretcher for compact storage
- Foot-end wheels fold flat for storage and snap back into place
- Fits squad bench post and wheel cups, for use as a second stretcher
- Includes two 5' quick release patient restraints

0107714 Model 107, Burgundy
0101070 Model 107, Orange

Model 107 features folding posts and wheels for more versatility



	Length	Height	Chair Height	Width	Depth	Storage Length	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	73	9	53	20	25	39	21	350	24	41x21x6
Metric	185	22	135	51	64	99	10	159	11	104x53x15

Model 107-B4

- Fixed head-end posts are support legs
- Fix posts are in position when chair is unfolded
- Fixed foot-end wheels allow easy maneuvering
- Folding foot-end handles provide excellent grip
- Can be used with Model 385-F carrying case
- Includes two 5' quick release patient restraints

0107715 Model 107-B4, Burgundy
0101074 Model 107-B4, Orange
0311762 Storage Case, Lettered, Burgundy
0311760 Storage Case, Lettered, Orange

Model 107-B4 has fixed posts. Fixed wheels, and folding handles at foot-end.



	Length	Height	Chair Height	Width	Depth	Storage Length	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	74	6	54	20	26	39	19	350	22	41x21x6
Metric	188	15	13	51	65	99	9	159	10	104x53x15

Model 107-C

- Fixed foot-end posts provide stability when the stretcher is used flat
- Fixed wheels at both ends permit it to be rolled like a chair or loaded into a specialty ambulance
- Folding handles at both ends enable it to be used as a standard pole stretcher
- Includes two 5' quick release patient restraints

0107716 Model 107-C, Burgundy
0101072 Model 107-C, Orange
0311762 Storage Case, Lettered, Burgundy
0311760 Storage Case, Lettered, Orange

Model 107-C features an adjustable backrest, wheels, and handles at both ends



	Length	Height	Chair Height	Width	Depth	Storage Length	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	73	6	54	20	26	39	25	350	30	41x21x6
Metric	185	15	13	51	65	99	11	159	14	104x53x15

Ferno Pole Stretchers

- Heavy duty aluminum poles for easy handling
- Heavy duty vinyl-coated nylon covers
- Nonslip handgrips

Model 108 Single Fold Pole Stretcher

- Simplest pole stretcher available
- Folds at full length

0101080 Model 108, Orange

0107735 Model 108, Burgundy



CE

	Open Length	Folded Length	Open Width	Folded Width	Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	81	na	23¼	na	1½	8	350	9	85x5x3
Metric	206	na	60	na	4	4	159	5	216x13x8

Model 108-A Single-Fold Pole Stretcher

- Rugged spreader bars and legs that provide support
- Folds at full length for easy storage
- Available in two lengths and colors

0101081 Model 108-A, 81" Orange

0107736 Model 108-A, 81" Burgundy

0101082 Model 108-A, 90" Orange

0107737 Model 108-A, 90" Burgundy

0311580 Storage case for model 108-A-81, Orange, Lettered

0311610 Storage case for model 108-A-90, Orange, Unlettered

Model 108-A



CE

	Open Length	Folded Length	Open Width	Folded Width	Height	Weight 81"/90"	Load Limit	Shipping Weight 81" / 90"	Shipping L,W,H
Imperial	81 / 90	na	21½	9¾	5½	12 / 13	350	9 / 15	93x7x5
Metric	206 / 229	na	55	25	14	5½ / 6	159	5 / 7	236x18x13

Model 108-AF Double-Fold Pole Stretcher

- Folds lengthwise and crosswise for compact storage
- Rugged spreader bars and legs provide support and comfort
- Available in both 81" (206 cm) and 90" (229 cm) lengths

0107982 Model 108-AF, 81" Orange

0107983 Model 108-AF, 81" Burgundy

0101084 Model 108-AF, 90" Orange

0107739 Model 108-AF, 90" Burgundy

0311470 Storage case for 108-AF-81, Lettered, Orange

0311490 Storage case for 108-AF-81, Unlettered, Orange

Model 108-AF can fold lengthwise and crosswise for compact storage



CE

	Open Length	Folded Length	Open Width	Folded Width	Height	Weight 81"/90"	Load Limit	Shipping Weight 81" / 90"	Shipping L,W,H
Imperial	81 / 90	40 / 45	21½	9¾	5¾	14 / 15	350	14 / 15	47x8x5
Metric	206 / 229	103 / 113	55	25	15	7 / 7	159	7 / 7	119x20x13

Model 108-F Double Fold Pole Stretcher

- Folds lengthwise and crosswise for compact storage

0107741 Model 108-F, Burgundy

0101086 Model 108-F, Orange

Model 108-F



CE

	Open Length	Folded Length	Open Width	Folded Width	Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H
Imperial	84	42	22	9¾	2	8	350	10	44x6x4
Metric	212	106	56	25	4	4	159	5	112x15x10