



Ferno Model 71 Basket Stretcher

- One piece basket stretcher
- High-density polyethylene shell
- Heavy-duty aluminum frame
- Molded-in runners and handholds
- Chemical, UV, rust, and corrosion resistant
- Brass grommets for quick use with Model 418-1 Adjustable Lifting Bridle
- Closed-cell foam pad is replaceable
- X-ray translucent for arms, legs, and chest
- Includes four restraints



0107100 Model 71, Basket Stretcher, One-Piece
0311800 Model 387-A, orange storage case for Model 71



Ferno Model 71-M Mining Non-Metallic Basket Stretcher

One piece mining version for use in subways and mining areas where electrical hazards or sparks could hamper rescue efforts. Plastic coated wood frame, plastic rivets, wooden runners, and hook-and-loop restraints.



0107104 Model 71-M, Mining Non-Metallic One-Piece Basket Stretcher

	Length	Width Length	Height	Weight	Load Limit	Shipping Weight 71	Shipping Weight 71M	Shipping L,W,H
Imperial	86	24	8	23	600	36	40	87x25x8
Metric	218	61	19	10	272	16	18	221x64x20

Ferno Model 71-S Split Basket Stretcher

- Two piece backpackable basket stretcher
- Has the same features as the one-piece Model 71 basket stretcher
- Separates in half for backpacking to a rescue site or compact storage

0107101 Model 71-S, Split Basket Stretcher, Two Piece
0311820 Model 387-B, orange storage case for model 71-S



	Length Length	Width Length	Height	Weight	Load Limit	Shipping Weight	Shipping L,W,H	
Imperial	86	24	8	23	600	38	47x27x11	
Metric	218	61	19	10	272	17	119x69x28	

Flotation Collar for Model 71 series



P/N 0313943
 • 45 lb. of bouyancy

Ferno Model 418-1 Lifting Bridle, Adjustable



0312390

Model 513-4 IV Pole



0085480



Ferno Model 108-D Disaster Litter

- The Model 108-D Disaster Litter is an economical, lightweight litter designed for disaster use
- Ideal for use by EMS services, fire departments, law enforcement, hospitals or other agencies that need a sturdy, compact and economical litter for mass casualty response
- Sturdy Ferno aluminum frame
- Easy-to-decontaminate polypropylene mesh cover
- Includes two attached polypropylene restraints
- Folds in half for storage

0107969 Model 108-D, Disaster Litter



CE

	Open Length	Folded Length	Open Width	Folded Width	Open Height	Folded Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	81	40	23	9	6	4	12	350	15	47x8x5
Metric	206	102	58½	23	15¼	10	5	159	7	119x53x13

Ferno Model 108-T Tactical Litter

- Designed for tactical, HazMat-CBRN Decon, USAR and disaster use
- Heavy duty, light weight poles, hinges, legs and spreader bars
- Sanoprene rubber end plugs keep contaminants from penetrating the poles
- Cover is durable, flame retardant, nonslip, chemical and mildew resistant polypropylene mesh
- Fluids flow freely away from patients during the decontamination process
- Easy-to-decontaminate polypropylene mesh cover
- Includes two attached polypropylene restraints
- Follows US military Spec MIL-L-49511B

0107970 Model 108-T, Tactical Litter, Olive Drab Green only



CE

	Open Length	Folded Length	Open Width	Folded Width	Open Height	Folded Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	92	46½	23	8	6	7	17	1200	21	49x12x8
Metric	234	118	58½	20	15¼	18	8	544	10	124x30x20

Ferno Model 135 Combination Stretcher/Body Bag

- Folds compactly and completely conceals the body
- Six durable, permanently blocked, formed handles
- Water-repellent top that closes with a sturdy slide fastener
- Two-ply, vinyl-coated nylon construction
- Resists stains and does not absorb moisture or bodily fluids
- Flame resistant
- Available in burgundy only

0310180 Model 135, Combination Flexible Stretcher/Body Bag, Burgundy



CE

	Length	Width	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	78	28	15½	350	22	91x6x5
Metric	198	72	7	159	10	231x15x13



Ferno Model 131 Flexible Stretchers

- Assists in moving patients through narrow openings or down steep staircases, and in handling unruly patients
- Solid-oak slat supports encased-in heavyduty, vinyl coated nylon
- Stretcher seams are electronically welded so blood and bodily fluids can't penetrate the surface
- Six carrying handles enable EMTs to maneuver patient
- Available model with built-in restraints that include chest, hip and individual leg straps
- Flame resistant
- Available in burgundy and orange



0310151 Model 131, Flexible Stretcher, Orange

0310161 Model 131, Flexible Stretcher, Burgundy



	<i>Length</i>	<i>Width</i>	<i>Weight</i>	<i>Load Capacity</i>	<i>Shipping Weight 131</i>	<i>Shipping Weight 137</i>	<i>Shipping L,W,H</i>
<i>Imperial</i>	78	28	15½	350	20	22	91x6x5
<i>Metric</i>	198	72	7	159	9	10	231x15x13

Ferno Model 137 Flexible Stretchers

- Same features as Model 131 plus built in restraints

0310200 Model 137, with Built-in Restraints, Orange

0310210 Model 137, with Built-in Restraints, Burgundy



Specifications the same as Model 131

Flexible Stretcher with built-in restraints >>



Ferno Model 132 Manta® Rescue Aid

- Compact and easy-to-store. Can be used to transfer bariatric patients, evacuate patients or cover patients during window extractions
- Ideal for rescue from confined areas
- Pockets at both ends hold backboard, patient's feet and head
- Fire retardant with holes at each end to allow for water drainage
- Available in burgundy, orange and olive drab green
- Folds like a blanket for easy storage



0382888 Model 132, Manta, Orange

0382969 Model 132, Manta, Burgundy

0383034 Model 132, Manta, Olive Drab Green

	<i>Length</i>	<i>Width</i>	<i>Weight</i>	<i>Load Capacity</i>	<i>Shipping Weight</i>	<i>Shipping L,W,H</i>
<i>Imperial</i>	78	45	4	800	7	17x11x7
<i>Metric</i>	198	115	2	363	3	43x28x18