

Millennia™ Boards

- The Millennium Board is made of a lightweight foam-filled construction
- All Millennium Backboards are X-ray translucent
- Can be used for water rescue
- Slick surface for easy transfer
- Minimal deflection for improved patient stabilization
- Wider, longer handholds for increased EMT control
- Compact tapered design for easy maneuverability in close confinement and storage

- Sealed construction for resistance to fluids, chemicals and temperature
- Ten X-ray translucent, polycarbonate speed-clip pins arranged for various strapping configurations
- Special mounting slots for the Ferno Model 445 Head Immobilizer
- Available in 16" (41 cm) and 18" (46 cm) widths
- Colors include burgundy, yellow, white, orange, and blue

18" Millennium Backboard

0107915 White
0107916 Blue
0107917 Orange
0107918 Yellow
0107928 Burgundy

16" Millennium Backboard

0107927 Burgundy
0107914 Yellow
0107913 Orange
0107911 Blue



Limited Lifetime Warranty

16" Board



	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	16	2	11	350	14	74x19x3
Metric	183	41	5	5	159	6	189x48x8

18" Board



	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	18	2	12	350	17	74x19x3
Metric	183	46	5	6	159	8	189x48x8



NAJO® Backboards

- **100% X-ray Translucent**
- **Impervious to Bodily Fluids and Bloodborne Pathogens**
- **All boards meet OSHA recommendations and are easy to decontaminate**
- **Lightweight and Buoyant**
- **Can float a patient in the water and can be used as a flotation device**
- **Guaranteed for Life**
ALL backboards have a limited lifetime warranty under normal usage.
- **One-piece, high-density, polyethylene plastic construction, with polyurethane foam fill**



Standard Colors

Najo Backboards are available in standard colors: yellow, orange, green, red, and blue.

Custom colors can be accommodated for larger orders (a minimum of 50 required).



Split Colors

Our boards are also available in a “split color” version, taking two of the colors shown above and applying them to opposite ends of the backboard.



Custom Graphics

- All NAJO Backboards (except the Disaster Board and NAJO Lite) can be ordered with custom graphics
- Place your name or logo on the board to identify it from other boards
- Graphics are permanent and cannot be removed
- Custom graphics have an extra cost



NAJO RediHold®



- Our top-of-the-line backboard
- Two handles at tapered end to allow for more stable transport
- Head area is smooth to insure all head immobilization devices adhere to the board
- Angled edges for log-rolling
- Able to float a weight of 155 lbs / 70 kg

16"

Top of the line



Limited Lifetime Warranty



	Yellow	Orange	Green	Red	Blue	White
No Pin	275101105	275101108	275101106	275101104	275101103	275101101
10 Pin	275101205	275101208	275101206	275101204	275101203	275101201
18 Pin	275101305	275101308	275101306	275101304	275101303	275101301

	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	16	1¾	16	600	20	73x18x2
Metric	183	41	4½	7	272	9	185x46x5

NAJO RediWide®

- 18" wide for rescuers needing "more backboard"
- 1 3/4" thick with angled edges for smooth log-rolls and painless patient placement
- Rectangular shape to allow patient placement at either end

18"



Limited Lifetime Warranty



	Yellow	Orange	Green	Red	Blue	White
No Pin	275301105	275301108	275301106	275301104	275301103	275301101
10 Pin	275301205	275301208	275301206	275301204	275301203	275301201
18 Pin	275301305	275301308	275301306	275301304	275301303	275301301

	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	18	1¾	16	600	20	73x18x2
Metric	183	46	4½	7	272	9	185x46x5

NAJO Lite Board

- Newest and most cost-effective X-ray translucent backboard
- Two handholds at the tapered end allow more stable transport
- Angled edges for log-rolling
- Central slot in the lower half of board to immobilize each leg

16"

Most cost effective



Limited Lifetime Warranty



	Yellow	Orange	Green	Red	Blue	White
No Pin	275201105	275201108	275201106	275201104	275201103	275201101
10 Pin	275201205	275201208	275201206	275201204	275201203	na
18 Pin	275201305	275201308	275201306	275201304	275201303	275201301

	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	16	1¾	14½	450	20	73x18x2
Metric	183	41	4½	7	204	9	185x46x5



NAJO RediBoard®

- Perfect water rescue board floats 275 lb / 125 kg
- Does not have a designated front, back, top or bottom, so patient can be positioned and immobilized faster
- 100 % X-ray translucent

16"



Water Rescue

Limited Lifetime Warranty



	Yellow	Orange	Green	Red	Blue	White
No Pin	275401105	275401108	275401106	275401104	275401103	275401101
10 Pin	275401205	275401208	275401206	275401204	275401203	275401201
18 Pin	275401305	275401308	275401306	275401304	275401303	275401301


	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	16	2¼	16	600	20	73x18x2
Metric	183	41	6	7	272	9	185x46x5

NAJO Disaster Board®

- Same great features as the RediBoard
- Sold only in sets of 35
- Only in yellow with no graphics and no pins
- Economic way to have backboards on hand when needed for a mass casualty disaster
- Does not have a designated front, back, top or bottom, so patient can be positioned and immobilized faster

ØDISASTERBOARD Set of 35, no pins, Yellow only

16"



Limited Lifetime Warranty



	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H	• Shipping information reflects one board
Imperial	72	16	2¼	15	600	20	73x18x2	
Metric	183	41	6	7	272	9	185x46x5	

NAJO Mega 1000™

- 1000 lb / 454 kg load capacity
- X-ray translucent backboard
- Features two tapered ends, allowing for extreme flexibility in transport

18"



Limited Lifetime Warranty



	Yellow	Orange	Green	Red	Blue	White
No Pin	275901105	275901108	275901106	275901104	275901103	275901101
10 Pin	275901205	275901208	275901206	275901204	275901203	275901201
18 Pin	275901305	275901308	275901306	275901304	275901303	275901301

	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	18	2¼	20	1000	20	73x18x2
Metric	183	46	6	9	454	9	185x46x5

NAJO Sports Board

- Specially designed for the professional athlete
- Large size allows athletes to be immobilized with equipment still on
- 100% X-ray translucent



20"



Limited Lifetime Warranty



	Yellow	Orange	Green	Red	Blue	White
No Pin	275801105	275801108	275801106	275801104	275801103	275801101
10 Pin	275801205	275801208	275801206	275801204	275801203	275801201
18 Pin	275801305	275801308	275801306	275801304	275801303	275801301

	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	80	20	2¼	27	1000	35	83x23x3
Metric	203	51	6	12	454	16	211x58x8

NAJO Water Rescue Backboard

- Specially designed for water and land rescue
- Lightweight and buoyant
- Tapered end without foam fill allows more effective submersion under patient in water rescue
- Design allows for easy placement during water rescue

275W01104 Water Rescue, No Pins, Red only

16"



Limited Lifetime Warranty
Available in Red Only



	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	16	2	14	450	20	73x18x2
Metric	183	40½	5.8	6.3	204	9	185x46x5

AquaBoard

- Designed and developed for the rescue of patients with suspected spinal injuries from swimming pools
- Also compatible for salt water rescue
- Neoprene Back Support attaches to the board surface and prevents the casualty from slipping once strapped to the board
- Instructions for patient positioning printed on the neoprene support
- Board floats up and under the casualty
- 100% X-ray translucent
- One year warranty

IE10309999 AquaBoard

17"



	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping H,L,W
Imperial	72	17	2¾	17½	400	23	78x17x4
Metric	183	43	7	8	182	10	201x43x10



Ferno Model 740 and 750 Series Phenolic Wooden Backboards

16"

- Constructed of 5/8" finished birch plywood
- Permanently coated with a phenolic resin to resist bodily fluids
- Optional runners make transport over rough terrain easy
- Available with or without speed-clip pins
- Unique Ferno EdgeGuard protects perimeter of backboard from splintering and cracking
- Reflective lettering on Ferno EdgeGuard increases nighttime visibility for extra safety
- Tapered design follows body contours for controlled handling in extrications or confined-space rescues
- Standard 16" width allows straps to be anchored close to the body, minimizing lateral movement during transport
- Fits inside standard wire or plastic basket stretchers
- Available in dark brown or tan
- One year warranty



Dark Brown Phenolic Back Board

0107745 Model 750, Plain

0107746 Model 751, with runners

0107747 Model 752, with pins

0107748 Model 753, with pins and runners

Tan Phenolic Back Board

0107722 Model 740, Plain

0107723 Model 741, with runners

0107724 Model 742, with pins

0107725 Model 743, with pins and runners



	Length	Width	Height w/runners	Height wo/runners	Weight w/runners	Weight w/o runners	Load Capacity	Shipping Weight	Shipping L,W,H
Imperial	72	16	1-5/8	1	18	15	350	21	75x16x3
Metric	183	41	4	2	8	7	159	10	191x41x8

Ferno Model 77 Short Wooden Backboard

13"

- Narrow board design for use in small autos or confined areas
- Multiple strap openings allow various strapping methods
- Notched head-end section enables head straps to be positioned according to patient size
- Constructed of 7-ply, 1/2" marine-grade plywood and finished with clear varnish on both sides
- Includes short wooden backboard, three color-coded 2" wide nylon patient restraints with automotive-type quick-release buckles, Adjusta-Pad neck roll and two head straps
- One year warranty



0107700 Model 77, Short Wooden Backboard Kit




	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping H,L,W
Imperial	32	13	1	4	350	9	34x14x4
Metric	81	32	2	2	159	4	86x36x10

Ferno Model 60 Aluminum Backboards



Durable, textured aluminum won't delaminate or chip and is easy to keep clean. Sturdy sheet aluminum is permanently fastened to a tubular aluminum frame forming a lightweight backboard.

- Multiple strap and handhold openings enable easy handling and various restraint configurations
- Folds in half to store in a variety of ambulance compartments
- Includes three 9' nylon patient restraints with quick-release metal buckles
- Available with three full-width runners mounted permanently to the underside of the backboard, enabling convenient grip and handling of the board
- Additional strap openings allow for greater flexibility in restraint configurations or for using cravats to immobilize the patient



CE

0106000 Model 60, with restraints

0106020 Model 60-2, with runners and restraints

0106040 Model 60-4, with runners and restraints

Model 60/60-2

	Open Length	Folded Length	Width	Open Height 60	Height with Runners 60-2	Folded Height 60	Height with Runners 60-2	Weight	Weight with Runners	Load Limit 60	Load Limit 60-2
Imperial	72	36	18	1	2¼	3	2½	12	15	350	350
Metric	183	91	46	3	5¾	7	6	5	7	159	159

Model 60-4

	Open Length	Folded Length	Width	Open Height	Folded Height	Height with Runners	Weight	Load Limit
Imperial	78¾	36	18	2¼	2½	2½	15	350
Metric	199	99	46	5¾	6	6	7	159

Ferno Model 61 Short Aluminum Backboard Kit

Narrow board design permits EMTs to apply short backboards in small autos or confined areas. Multiple strap openings allow various strapping methods depending on the situation and the condition of the patient.

- Notched head-end section enables head straps to be positioned according to patient size
- Constructed of sheet aluminum fastened to an all-aluminum frame
- Includes short aluminum backboard, three color-coded 2" wide nylon patient restraints with automotive-type quick-release buckles, Adjusta-Pad neck roll and two head straps



CE

0106100 Short Aluminum Backboard Kit

	Open Length	Width	Open Height	Weight	Load Limit	Load Limit
Imperial	32	13	1	3	350	500
Metric	81	32	2	2	159	227