



## Model 65 Scoop™ Stretcher

- Can uncouple at one or both ends of the stretcher to gently scoop up the patient using a scissor-type closing motion
- Ideal for hip injuries and accident victims found on the street, because it supports a patient in the position found, reducing risk of further injury
- Narrow foot-end frame for handling in confined areas
- Telescoping tubes allow length adjustments to fit patients of various heights
- Open-center design permits patient to be X-rayed
- Folds in half at foot end for compact storage
- Can be used with Model 385-C Storage/Carrying Case (#0311700)
- **Model 65-PIN Scoop™ Stretcher** includes eight speed-clip pins arranged for various strapping configurations (#0000325)



Telescoping poles  
for precise fit

- 0000384** Model 65, with Burgundy Restraints and metal buckles  
**0000385** Model 65, with Orange Restraints and metal buckles  
**0000365** Model 65, with Orange polypropylene straps and plastic buckles  
**0000386** Model 65, with Burgundy polypropylene straps and plastic buckles  
**0000325** Model 65, with Black polypropylene straps and speedclips



	Open Length	Folded Length	Width	Open Depth	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping L,W,H
<b>Imperial</b>	65	47	17	3	4	16	350	26	67x18x3
<b>Metric</b>	166	120	43	7	9	7	159	2	170x46x8

## Ferno Model 65 EXL Scoop Stretcher

Designed to uncouple either end or both ends to gently slide scoop under victim to support patient in the position found with minimal disturbance to reduce risk of further injury


- Thermally-treated polymer construction does not get too hot or cold, allowing for greater patient comfort
- Oval profile and angled handles provide ergonomic ease of use
- Twin safety lock for smooth locking and unlocking
- Head section recessed for proper cervical alignment
- Four length options for patients of different heights 65", 70", 74", 79"
- Narrow foot-end frame for handling in confined areas
- Easy-to clean surface is impervious to fluids
- Folds to 47" for compact storage

### Shipped from Italy

- 0107957** 65 EXL with pins, three #430 restraints, 5' two piece metal buckles  
**0107958** 65 EXL with pins, three #430 restraints, 7' one piece metal buckles  
**0107959** 65 EXL with pins, three #430P restraints, 5' two piece plastic buckles  
**0107960** 65 EXL with pins, three #430P restraints, 7' one piece plastic buckles

### Shipped from Wilmington

- 00065EXL** 65 EXL with pins, three #430P restraints, 5' two piece plastic buckles  
**0107991** 65 EXL Scoop with pins and without restraints



Includes restraints  
(burgundy)

X-Ray translucent

Non-Magnetic

Twin Lock



	Open Length	Folded Length	Width	Open Depth	Folded Depth	Weight	Load Limit	Shipping Weight	Shipping L,W,H
<b>Imperial</b>	65	47	17	3	4	18	350	22	67x18x3
<b>Metric</b>	165	120	43	7	8	8	159	10	170x46x8