

# Ferno Model 770 Fastrap™ Quick Restraint System for Backboards and Scoops



- Ten-point, six-strap system used to secure patients to backboards
- Nylon with hook-and-loop closures
- Size adjustable from pediatric to adult
- Compatible with the Millennia Backboard, Phenolic Wooden Backboard Series, 65 Scoop, and 65 EXL Scoop Stretchers

**0313916** Fastrap Quick Restraint System



CE

## Ferno Model 430 Cot Restraint Straps

- One- or two-piece, primary-use patient restraints made of 2" nylon webbing
- Quick-release, automotive-type, chrome-plated metal buckles
- Rust and corrosion resistant
- Push-button operation allows secure fastening and easy release
- Hemmed tab prevents strap from being pulled free of buckle



### Two-Piece

	Black	Burgundy	Orange	Gray
<b>5'</b>	0313892	0313890	0313891	0313893
<b>7'</b>	0313911	0313913	0313912	0312620
<b>9'</b>	0313305	0313306	0312560	na

### One-Piece

	Black	Burgundy	Orange	Gray
<b>7'</b>	0313910	0313888	0313887	0312610
<b>9'</b>	0313899	0313909	0313889	na

CE

## Ferno Model 430P Restraint Straps

- One- or two-piece, secondary-use patient restraints made of 2" polypropylene webbing
- Side-release, snap-lock, plastic buckle tolerates aquatic situations without water damage
- Squeezable outer prongs allow secure fastening and easy release
- Hemmed tab prevents strap from being pulled free of buckle



### Two-Piece

	Black	Burgundy	Orange	Gray
<b>5'</b>	0313797	0313918	0313917	0313796
<b>7'</b>	na	0313800	0313919	0313795
<b>9'</b>	0313805	0313806	na	na

### One-Piece

	Black	Burgundy	Orange	Gray
<b>7'</b>	0313894	0313895	0314011	0313794
<b>9'</b>	0313790	0313793	0313789	na

CE



## Ferno Model 436 and 436P Swivel Speed-Clip Restraints

Two-piece patient restraints made of 2" nylon/polyester webbing with quick-release, automotive-type, coated metal buckle that resists rust and corrosion (Model 436), or polypropylene webbing with plastic buckle (Model 436P)

- Hemmed tab prevents strap from being pulled free of buckle
- Swivel speed-clips for quick and easy attachment to backboard with pins
- 5' long, 2 piece, available in burgundy and orange

### Model 436



### Model 436P



### Two-Piece Only

- 0313079** Model 436 with Metal Buckle, Set of 3, Orange  
**0313082** Model 436 with Metal Buckle, Set of 3, Burgundy

- 0313809** Model 436P with Plastic Buckle, Set of 3, Burgundy  
**0313785** Model 436P with Plastic Buckle, Set of 3, Orange



## Model 439 and 439P Backboard Strap Kits with Bag

- Each kit includes a set of three straps and a storage bag
- One piece restraints with swivel speed clip
- 7' long only, available in black, burgundy, and orange
- Metal or Plastic buckles



### One-Piece Only

- 0313083** Model 439 Strap Kit, with metal buckles, Set of 3, Black  
**0313084** Model 439 Strap Kit, with metal buckles, Set of 3, Burgundy  
**0312660** Model 439 Strap Kit, with metal buckles, Set of 3, Orange

- 0313747** Model 439P, with plastic buckles, Set of 3, Burgundy  
**0313746** Model 439P, with plastic buckles, Set of 3, Orange

