



## Ferno Model 78 Pedi-Pac™

- A spine immobilizer and restraint system for children from 28" to 54" (71 to 137 cm) tall, weighing 20 to 90 lb (9 to 41 kg)
- Built-in fastening loops connect to existing cot straps for safety during transport. Individual leg restraints allow one leg to be immobilized while EMT attends to other leg
- Adjustable head support with ear openings for monitoring fluid drainage
- Sewn-in lifting handles at both ends for easy handling in confined areas
- Replaceable color-coded straps for easy identification
- Unit comes complete with head and chin straps and carrying case



**0313065** Model 78, Pedi-Pac



	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
<b>Imperial</b>	47	8	1	6	90	10	49x10x4
<b>Metric</b>	120	20	4	3	41	5	125x25x10

## Ferno Model 125 Kendrick Extrication Device, KED500

- Improved means of immobilizing and extricating patients from auto accidents or confined spaces
- The wrap-around design provides vertical rigidity for maximum support of the spine, neck, and head during extrication
- Immobilizes torso, head and neck minimizing risks of further injury
- X-rays or advanced life support procedures can be performed with the K.E.D.® in place
- Heavy duty, vinyl-coated nylon is durable and easy to clean
- Case is enlarged to accommodate extrication collars
- Built-in handles enable rescuers to get a firm grip on the patient and K.E.D.® as one unit
- Adapts to a hip and pelvic splint by simple inversion
- Adjustable fold-back sides permit easy access to patient's chest
- Can be used for multiple patient sizes, including pediatric patients, children, and pregnant women

**0313676** Model 125, KED500



	Length	Width	Height	Weight	Load Capacity	Shipping Weight	Shipping L,W,H
<b>Imperial</b>	33	32	6	8	500	10	34x10x5
<b>Metric</b>	83	80	14	3	227	5	86x25x13