

Ferno Model 445 Universal™ Head Immobilizer

Composed of two plastic-coated, closed-cell foam head supports, one universal attachment base and two durable head straps. Base plate fits all Ferno wooden, aluminum and plastic backboards, and scoop stretchers. Head supports can be used on standard backboards, or turned around for use on the Scoop Stretcher

0313855 Model 445, Universal Head Immobilizer

- Unique head strap is designed to conform to the patient's forehead
- Large ear holes for monitoring the patient's ear canal
- Minimal interference with X-ray, MRI, or CT scanning procedures
- Waterproof plastic coating makes it easy to clean and prevents bacterial growth within components
- Will not absorb blood or bodily fluids; easily sanitized between uses
- Meets FMVSS 302 flammability rating



Compatible with all Ferno backboards and scoops





Ferno Model 455 HeadHugger® Head Immobilizer

- Low-cost method for effective head immobilization
- Heavy-grade, water-resistant corrugated material
- Base plate secures to the board with adhesive tabs
- Independent support panels fold around the patient's head and are held in place with hook-and-loop fasteners
- Adhesive tape secures the patient's head to the Head-Hugger
- Does not interfere with CT scanners or MRIs
- Ideal for use in mass casualty accidents
- Easy access to monitor blood loss or fluid drainage
- Dual head supports can be molded to fit the contour of each patient
- Sold in sets of 10

0831919 Model 455 HeadHugger, sold in sets of 10



Ferno Model 449 WizLoc® Cervical Extrication Collar

- WizLoc provides a collar capable of fitting 90% of all patients
- Unique self-centering mandibular support allows custom fit to the patient
- Designed for use on patients 2 years old through adult
- On-patient sizing for greater accuracy
- Three position occipital support with keyhole lockpin
- Reusable up to 10 times
- Can be decontaminated
- High visibility colors for safety and identification
- X-ray translucent and MRI compatible
- Also available in military OD green (#0819759, Model 449-M)

0819750 WizLoc® Cervical extraction Collar, sold in sets of 3



Mandibular Support



Occipital Support





Philly® One-Piece Cervical Collar

- Very simple to apply
- Application in a seated position is quick and easy
- Lightweight, thermo-formed and pre-molded chin support design
- Enlarged trachea opening
- Color-coded sizing
- Radiolucent
- Available in tall neck (green), regular (orange), short neck (blue), low (purple), pedi (aqua) or infant (pink)

One-Piece Cervical Collar

IE11309901 Regular, Orange

IE11309908 Infant, Pink

IE11309909 Pediatric, Aqua

IE11309910 Low, Purple

IE11309911 Short Neck, Blue

IE11309912 Tall Neck, Green



Philly® Patriot Cervical Collar

- Adult collar adjusts to four different heights
- Pediatric collar adjusts to three different heights
- Injection-molded polypropylene
- Patented chin support
- Oversized trachea opening
- Open rear panel
- Essentially flat design
- Nasal cannula holders
- Radiolucent

IE11009907 Adult, White and Blue

IE11109906 Pediatric, White and Blue



Philly® Philadelphia Red Cervical Collar

- Conforms to the spine to support the cervical arch with a pre-formed back
- Supports the head carefully without tilting it from side to side and eliminates rotation with a head immobilizer attachment
- Available in six sizes: infant, child, short neck, adult, tall neck, and extra-large

IE11200306 Child, Red only

IE11200311 Short Neck, Red only

IE11200307 Adult, Red only

IE11200303 Extra Large, Red only

IE11200308 Infant, Red only

IE11200312 Tall Neck, Red only

