



Model 441, 443, 444 FernoTrac® Traction Splints

- Provides in-line mechanical traction
- Relieves pain and reduces the possibility of further vascular and nerve damage
- Simple, adjustable, wrap-around ankle hitch designed to allow EMTs to check pedal pulses
- Low-profile ischial pad permits splint to be applied with minimum movement of the patient's leg
- Heel stand folds flat for storage, and instantly locks into place for support of the injured leg

0819944 Model 444, Adult Traction Splint

0819945 Model 443, Pediatric Traction Splint

0819946 Model 441, Combination 443 plus 444

0313031 Model 441-1, Carrying Case Only



	<i>Min Length</i>	<i>Max Length</i>	<i>Width</i>	<i>Open Height</i>	<i>Folded Height</i>	<i>Weight</i>	<i>Shipping Weight</i>	<i>Shipping L,W,H</i>
Imperial	36	54	11	9	3	6	18	40x12x5
Metric	90	137	27	22	8	3	4	102x31x13

Sager® Traction Splints

- Shows the exact amount of traction being applied
- Bilateral Model fits within the adult body silhouette
- One splint provides traction splinting for single or bilateral fractures
- Permits the traction to decrease as the spasm releases
- One person application
- Nothing to assemble, "always ready to use"
- Compatible with anti-shock trousers
- Lightweight stainless construction
- No heel stands
- Incorporates an articulating ischial perineal cushion that is comfortable and conforms to the pelvis, simplifying pelvic procedures
- Unilateral and bilateral models fit a child of 4-years up to an adult over 7' tall (2.13 mts)
- Storage Bag and application instructions included

SP3010 Unilateral, for children of 4 years up to adults

SP3040 Bilateral, for children of 4 years up to adults

SP3000 Infant Bilateral, from infants up to children 6 years of age



	<i>Min Length</i>	<i>Max Length</i>	<i>Width</i>	<i>Open Height</i>	<i>Folded Height</i>	<i>Weight</i>	<i>Shipping Weight</i>	<i>Shipping L,W,H</i>
Imperial	32	47	11	8	3	6	12	38x21x4
Metric	82	120	27	19	8	3	5	97x53x10



Kendrick Traction Device (KTD)

- Eliminates the need to roll the patient or raise the leg during application
- The absence of an ischial bar eliminates problems of needless pressure
- Allows to work around hip or groin
- Small enough to be carried in a first aid kit
- Adjusts to adult and pediatric sizes
- Compatible with anti-shock trousers
- Available in orange and black

0819766 KTD, Orange

0819760 KTD, Black



Ferno Model 190 Vacuum Splint Extremity Kit

- Ideal for dislocated shoulders, suspected fractures of arms, legs, and joints or infant immobilization
- Kit includes three splints (small, medium, and large), compact vacuum pump, repair kit and carrying case.

0382885 Model 190 Complete Extremity Kit, Orange

Splints also available individually

0382910 Model 100 Short arm, Small, Orange

0382912 Model 120 Long arm, Medium, Orange

0382914 Model 140 Leg, Large, Orange

0831965 Replacement Vacuum Pump

0382916 Carry Case Only, for Model 190



CE

Ferno Model 192 Full-Body Vacuum Mattress

- Immobilizes the patient by conforming to any body shape
- Three patient restraints with buckles provide maximum patient security
- Six carry handles allow the vacuum mattress to be used as a stretcher
- Includes removable, contoured pillow and compact vacuum pump

0382436 Model 192, Full Body Vacuum Mattress, Orange

0831965 Replacement Vacuum Pump



CE