



## Kendrick Traction Device (KTD)

- Eliminates the need to roll the patient or raise the leg during application
- The absence of an ischial bar eliminates problems of needless pressure
- Allows to work around hip or groin
- Small enough to be carried in a first aid kit
- Adjusts to adult and pediatric sizes
- Compatible with anti-shock trousers
- Available in orange and black

**0819766** KTD, Orange

**0819760** KTD, Black



## Ferno Model 190 Vacuum Splint Extremity Kit

- Ideal for dislocated shoulders, suspected fractures of arms, legs, and joints or infant immobilization
- Kit includes three splints (small, medium, and large), compact vacuum pump, repair kit and carrying case.

**0382885** Model 190 Complete Extremity Kit, Orange

**Splints also available individually**

**0382910** Model 100 Short arm, Small, Orange

**0382912** Model 120 Long arm, Medium, Orange

**0382914** Model 140 Leg, Large, Orange

**0831965** Replacement Vacuum Pump

**0382916** Carry Case Only, for Model 190



CE

## Ferno Model 192 Full-Body Vacuum Mattress

- Immobilizes the patient by conforming to any body shape
- Three patient restraints with buckles provide maximum patient security
- Six carry handles allow the vacuum mattress to be used as a stretcher
- Includes removable, contoured pillow and compact vacuum pump

**0382436** Model 192, Full Body Vacuum Mattress, Orange

**0831965** Replacement Vacuum Pump



CE