



All ORO₂ Regulators offer:

- Strong uni-constructed body
- Low profile gauge
- CGA 870 pin index connection
- Two DISS high flow ports
- 0 - 3,000 psig Gauge
- Fixed orifice flow control with window-style flow reading
- Flow specifications range from 0-25 lpm
- Settings of ¼, ½, 1, 2, 3, 4, 6, 8, 10, 15 and 25
- Regulator provides a gas regulated pressure of 57 psig at 2000 psi
- Flow accuracy of plus or minus 10% of reading
- "C"-ring clamp aids in retaining the tee-handle
- Five year warranty and a voucher for a free spare replacement gauge, shipped with the unit when purchased

0819740 Model 816 - Aluminum Body, Brass Component
ORO₂ Oxygen Regulator, 0 - 25 lpm, 12-position flow control with barbed outlet and two DISS outlets with safety check valves, low profile, booted gauge.



Model 816

0819741 Model 817 - Aluminum Body, Brass Component
ORO₂ Oxygen Regulator, same as Model 816, except with international 0 - 4000 psi booted pressure gauge.



Model 817

0819742 Model 818 - Aluminum Body, Brass
ORO₂ Oxygen Regulator, 0 - 25 lpm, 12-position flow control with barbed outlet only, low profile, booted gauge.

0819743 Model 818d (not shown)

ORO₂ Oxygen Regulator, same as Model 818, except with a 12-position flow control with DISS outlet only.



Model 818

0819744 Model 819 - Aluminum Body, Brass Component
ORO₂ Oxygen Regulator, same as Model 818, except with international 0 to 4000 psi booted pressure gauge.



Model 819

0819567 Model 804 - Aluminum Body, Brass Component
ORO₂ Oxygen Regulator, CGA 540 DISS inlet for 'H' and 'M' size cylinders, one DISS 40 to 60 psi outlet with safety check valve UL listed 0 - 4000 psi pressure gauge with boot.



Model 804



FERNO₂ BRAVO₂™ Brass Oxygen Regulator

All BRAVO₂™ Regulators offer:

- Bodies constructed from a solid piece of brass
- A uni-body style regulator
- Internal brass components will not burn under 10,000 psi.
- Nickel-plated regulator
- Corrosion-resistant finish
- Available in CGA 870 and 540 configurations

0819745 Model 820

BRAVO₂™ Oxygen Regulator, 0 - 25 lpm, 12-position flow control with barbed outlet and two DISS outlets with safety check valves, low profile, booted gauge.



Model 820

0819746 Model 821

BRAVO₂™ Oxygen Regulator, same as Model 820, except with international 0 - 4000 psi booted pressure gauge.



Model 821

0819747 Model 822

BRAVO₂™ Oxygen Regulator, 0 - 25 lpm, 12-position flow control with barbed outlet only, low profile, booted gauge.



Model 822

0819748 Model 823

BRAVO₂™ Oxygen Regulator, same as Model 822, except with international 0 - 4000 psi booted pressure gauge.



Model 823