



## **FERNO<sub>2</sub> Disposable Bag Valve Mask Resuscitators (BVM)**

- Transparent re-breathing valve assembly allows for clear viewing of valve function
- Latex free
- Lightly-textured PVC bag for sure grip and comfort
- 360° patient valve swivel connection
- Adult size will easily collapse for compact storage
- Comes complete with mask reservoir bag and O<sub>2</sub> tubing
- Bag Valve Mask Resuscitator can be used in conjunction with FERNO<sub>2</sub> HFOC Model 301 to achieve 100% oxygen delivery
- 40 cm H<sub>2</sub>O pressure-relief valve standard on pediatric and infant models

**RB110005** Adult Bag Mask

**RB111009** Pediatric Bag Mask

**RB112008** Infant Bag Mask



## **FERNO<sub>2</sub> Silicone Resuscitators**

- Silicone material provides excellent bag re-expansion
- Superior feel of lung compliance helps reduce risk of over inflation
- Sized for adult, pediatric and infant patient ventilation (5.5-55 lb)
- Integrated intake/oxygen valve for ease of oxygen delivery
- Easy to disassemble/assemble for efficient cleaning
- See-through mask with unique multi-function mask cover
- 40 cm H<sub>2</sub>O pressure-relief valve standard on pediatric and infant models
- BVMs can be used in conjunction with FERNO<sub>2</sub> HFOC Model 301 to achieve 100% oxygen delivery
- Comes complete with mask reservoir bag and O<sub>2</sub> tubing

**RB121005** Adult Complete in carton

**RB122004** Child Complete in carton

**RB123008** Infant Complete in carton



## **FERNO<sub>2</sub> Reusable Silicone Face Masks**

- Contoured anatomically-shaped cushion
- Gas-tight low-pressure seal
- Hypo-allergenic materials
- Little or no blow-by during oxygen administration
- Available in pediatric and adult sizes

**RB132080** Reusable Face Mask, size 0, Infant

**RB1320081** Reusable Face Mask, size 1, Infant

**RB1325042** Reusable Face Mask, size 2, Infant/Child

**RB1330043** Reusable Face Mask, size 3, Child

**RB1330054** Reusable Face Mask, size 4, Adult

**RB1330055** Reusable Face Mask, size 5, Adult

