

- Reusable mask used for mouth to mask ventilation
- Ready to use right out of its rugged storage case
- Latex-free
- Comes complete with filtered inlet valve, O₂ inlet port and head strap

RO2820 FERNO₂ CPR Pocket Mask



CE

Ferno Model 672 CPR Board

- For CPR training and CPR administration by providing a firm surface
- Helps to properly position a patient to ensure a constant, open airway
- Head cup maintains proper alignment of patient's head
- Contoured sides make lifting the patient easier
- Measures 23¼x17x¼ (4½ lb)

0831820 CPR Board



FERNO₂ ST Fiber Optic Laryngoscopes

- Blades are reusable and can be sterilized
- ST Fiber Optic Blade eliminates the need to change or detach the bundle during the lifetime of the instrument

ST Fiber Optic Blades

0819594 Miller 1: 95 mmL, 13 mmW

0819595 Miller 2: 152 mmL, 13 mmW

0819596 Miller 4: 207 mmL, 15.9 mmW

0819597 Macintosh 2: 115 mmL, 13 mmW

0819598 Macintosh 3a: 145 mmL, 15 mmW

0819601 Spare Bulb: 2.5V, Halogen/Zeon

Laryngoscope Handle:

- ISO 2000 Connection
- Green System Compatible
- Stainless Steel Grip with Sanalon connection
- Halogen/Zeon bulb for a brighter high-intensity light
- The bulb is designed into the ISO 2000 handle for protection and eliminates the need to disassemble the blade before cleaning

0819599 Handle



Kit Model 425 includes:

- Carrying Case
- Miller 1, 2 & 4
- Mac 2 & 3a
- Handle
- Spare bulb (size "C" batteries not included)

0819593 Kit Model 425



Individual Laryngoscope Blades

Macintosh

- AM14040000** Macintosh Fiber Optic SS, size 0
- AM14040001** Macintosh Fiber Optic SS, size 1
- AM14040002** Macintosh Fiber Optic SS, size 2
- AM14040003** Macintosh Fiber Optic SS, size 3
- AM14040004** Macintosh Fiber Optic SS, size 4

Miller

- AM14050000** Miller Fiber Optic SS, size 0
- AM14050001** Miller Fiber Optic SS, size 1
- AM14050002** Miller Fiber Optic SS, size 2
- AM14050003** Miller Fiber Optic SS, size 3
- AM14050004** Miller Fiber Optic SS, size 4
- AM14010999** Handle Fiber Optic Laryngoscope, Pediatric
- AM14011199** Handle Fiber Optic Standard or Pediatric



Berman Airways

- Molded of semirigid, non-toxic, polyethylene
- Two-channel design allow for safe insertion or withdrawal of a suction catheter and ensures an unobstructed pathway for ventilation
- Kit includes six sizes: infant (43mm), child (60mm), small adult (80mm), medium adult (90mm), large adult (100mm), and extra-large adult (110mm)

AM0001 Complete Kit with Case, 6 per box



Guedel Emergency Airway Kit

- Combines six disposable, individually numbered Guedel Airways
- Plastic pouch for easy storage and access
- Large lumen offers free exchange of expired air and rescue breaths
- Each airway is contoured for maximum comfort
- Kit includes six sizes: infant (50mm), small child (60mm), child (70mm), small adult (80mm), medium adult (90mm) and adult (100mm)

AM21009999 Guedel Emergency Airway Kit, 6 per box





Color-Coded Guedel Airway Kit

- Six single use, color coded to identify correct size for rapid treatment of patients
- Features a rigid bite block to help prevent patient from closing off the airway
- Latex free

AM21019999 Complete Kit in a Box, set of 6

Complete set comes in plastic bag with standard sizes: infant (50mm), small child (60mm), child (70mm), small adult (80mm), medium adult (90mm) and adult (100mm).



Latex-Free Nasopharyngeal Airways

- Robertazzi style airway with fixed flange
- Made of neoprene and designed to minimize entry trauma and increase patient comfort
- For use where patient has a blockage of normal anatomical airway

AM1520 Latex-Free Naso Airway Kit of 9

AM15219999 Latex-Free Naso Airway Kit of 6

Kit of 9 includes nine sizes along with water-soluble sterile lubricant packs: 20, 22, 24, 26, 28, 30, 32, 34, and 36 fr. Also available individually. Kit of 6 excludes: 20, 34, and 36 fr.



FERNO₂ PosiTube Endotracheal Detection Device

- Ideal device to differentiate between the patient's trachea and esophagus
- Designed to verify proper placement of the endotracheal tube or double lumen airway

RO6200 PosiTube Endotracheal Tube Placement Check Device



FERNO₂ Grip-ET® Endotracheal Tube Holder

- Works with any ET tube from sizes 4.5 to 10mm
- The patented clip slides into place, locking quickly
- The wide open base leaves the mouth open for visualizing and suctioning the mouth
- The twist lock pivots hold the clip in for added security
- The Velcro® strap wraps around the head for a quick and easy hold

FA0009 Grip-ET® with Endotracheal Tube Holder Bite Block

