



Model 5 Intensive Care Transport

- The Model 5 is designed for transport of critical patients in Bell 212/412 Helicopters
- Unit is tested and certified to comply with FAA TSO 39b
- Can be installed without aircraft modifications by utilizing base adapters attached to standard floor fittings
- Installation requires less than 10 minutes in most aircraft



Model AV-5 Model 5, Intensive Care Transport System

Contact Ferno Aviation, Inc. for more details.

Please know the make, model, and type of aircraft that will be used for transport.

Each Unit is Equipped With:

- Ferno Model 12-2 Air Transport Stretcher with quick-release
- Shoulder harness attached to the unit's frame structure
- Storage drawers and compartments for drugs and disposables
- Fold-out instrument tray
- Cargo rings for easy tie-down
- Base interlock activated by fold-out carrying handles
- Pull-out swivel shelf and quick release mounts for one each LifePac 10 Cardiac Monitor/Defibrillator (LifePac 10 sold separately)
- Storage compartment with cylinder bottle lock for four Luxfer "D" size aluminum oxygen cylinders (oxygen bottles and regulators sold separately)
- Attachment point for one each LSP 3000 Automatic Ventilator (LSP 3000 sold separately)
- Attachment point for one each Matrix Emergency Suction unit with case (Matrix Emergency Suction unit sold separately)
- Shelf enclosure for one each Thomas Advanced Life Support Pack (Thomas Advanced Life Support Pack sold separately)
- Transporter (shown in picture) for the Model 5 is only for transport of the Model 5 to and from the aircraft. The Transporter is not intended for transport of the Model 5 with a patient

	<i>Length</i>	<i>Width</i>	<i>Height</i>	<i>Structural Framework and Stretcher</i>	<i>Medical</i>	<i>Miscellaneous Storage</i>
<i>Imperial</i>	73.87	20½	16¼	140	85	115
<i>Metric</i>	187.62	57.70	41.27	63.63	38.63	52.27

Model 46 Passenger Loading Chair

- Designed to transport a person through narrow passageways, including most commercial aircraft aisles, and other areas where space is limited
- Has a sturdy lightweight aluminum frame designed for easy loading by airline attendants
- Has a heavy duty vinyl coated nylon cover and seat pad
- The passenger loading chair comes with two upper body restraints and one leg restraint

Model AV-46 Model 46, Passenger Loading Chair 176

Contact Ferno Aviation, Inc. for more details.

Please know the make, model, and type of aircraft that will be used for transport.



	<i>Length</i>	<i>Width</i>	<i>Height</i>	<i>Weight</i>
<i>Imperial</i>	33¾	14½	52¼	17
<i>Metric</i>	85	37	133	7.71



The Columbia Model 8020-R RescueMate™

- Patient loading chair for use by fire and EMS personnel in both aircraft and ground vehicles
- Manufactured with ABS plastic, which meets FAR fire requirements
- Four-wheel brake system
- Telescoping lift handle
- Color-coded patient restraints
- Available in Yellow and Gray

AV-8020-R Model 8020-R, RescueMate, Yellow

Contact Ferno Aviation, Inc. for more details.



The Columbia Model 8020-T TravelAide™

- The same features as the Model 8020-R but in dark gray
- Both Columbia Chairs can be used as a passenger loading chair in airports, on buses, trains, leisure boats, cruise ships, and recreational vehicles

AV-8020-T Model 8020-T, TravelAide, Gray

Contact Ferno Aviation, Inc. for more details.

	<i>Extended Length</i>	<i>Folded Length</i>	<i>Extended Width</i>	<i>Folded Width</i>	<i>Extended Height</i>	<i>Folded Height</i>	<i>Weight</i>	<i>Load Capacity</i>
Imperial	27	12½	16	16	42	42	36	300
Metric	69	32	41	41	107	107	16.3	136.1

Model 8010 AisleMaster™ Boarding Chair

- Solid welded frame
- Flip-up armrests for easy transfer of passenger/patient
- Color-coded security belts
- "Deadman" brake handle
- Telescoping lift handles
- Adjustable headrest

Contact Ferno Aviation, Inc. for more details.

AV-8010 Model 8010, AisleMaster Boarding Chair



	<i>Extended Length</i>	<i>Folded Length</i>	<i>Extended Width</i>	<i>Folded Width</i>	<i>Extended Height</i>	<i>Folded Height</i>	<i>Weight</i>	<i>Load Capacity</i>
Imperial	27	na	12	na	42	na	33	na
Metric	69	na	30	na	107	na	15	na

Model 8030 TransportMate™ Boarding Chair

- Manufactured with ABS plastic, which meets FAR fire requirements
- Has folding push handles and ball bearing casters for easy maneuverability
- Four-wheel brake system
- Positive locking mechanism to prevent accidental folding
- Security belts

AV-8030 Model 8030, TransportMate Boarding

Contact Ferno Aviation, Inc. for more details.



	<i>Extended Length</i>	<i>Folded Length</i>	<i>Extended Width</i>	<i>Folded Width</i>	<i>Extended Height</i>	<i>Folded Height</i>	<i>Weight</i>	<i>Load Capacity</i>
Imperial	27.4	33	16	16	32	7	16	na
Metric	70	84	41	41	82	18	7.3	na