



Orion Bacteriostatic Blanket

- 100% liquid proof welded blanket
- Warm, comfortable, and light
- Washable at 90° C
- 110 x 190 cm
- Fabric-dark blue jersey, bacteriostatic (liquid proof on both sides)
- Fireproof (all materials classified M1)
- Can be sanitized with disinfectant



Vaccum Splints

- Soft and flexible heavy-duty material
- Removable straps for full and easy disinfection
- X-Ray translucent
- Color orange / white grey
- Max vacuum pressure - 800 hPa



Vacuum Mattress 94 ART



- PVC coated polyester sheet 750 gr/M2
- 3 (50 mm) restraints with buckles to provide maximum patient security
- 8 carry handles
- Groundsheet an option
- Colors: Blue and Orange
- EN 1865



CE

Vacuum Mattress 99 ART

- PVC coated polyester sheet 750 gr/M2
- Anatomical shape
- 3 (50 mm) restraints with quick release buckles to provide maximum patient security
- 8 ergonomic carry handles
- Groundsheet an option
- Temperature range 20° - 80° C
- Colors: Blue and Orange
- EN 1865



CE

Vacuum Mattress 99 ART PP

- PVC coated polyester sheet 750 gr/M2
- Total anti-reaction
- Anatomical shape
- Quarter turn valve
- 100% X-Ray translucent
- 3 (50mm) restraints with quick release buckles to provide maximum patient security
- 12 ergonomic handles for easy carrying
- 1 more handle at both ends
- Temperature range 20° - 80° C
- Groundsheet included
- Colors: Blue and Orange



CE