



Utila Chair Model ALS103

Our standard chair model is also the most popular one among the customers. All Utila chairs are delivered with one lap restraint. The attached picture shows the Utila chair ALS 103 with a shoulder harness restraint, which is an option.

- Strong aluminum construction
- White powder coated finish
- Anatomically shaped seat and backrest
- Washable synthetic leather, can be disinfected
- 2 rear telescoping handles
- 2 rear folding handles
- 2 front telescoping handles
- Folding armrests
- Telescoping footrest
- 2 stationary and two swivel wheels with brakes
- Equipped with 2 pins for locking into floor mount Utila model BH 501



	<i>Length Handles Extended</i>	<i>Length Handles Retracted</i>	<i>Width</i>	<i>Height</i>	<i>Seat Height</i>	<i>Weight</i>	<i>Load Limit</i>
<i>Metric</i>	150	80	57	110	50	23.5	150

Utila Chair Model ALS300

This chair is similar to our standard model but weighs only 15 kg. Therefore it is easy to handle though it is strong and stable.

- Standard features:
- Strong aluminum construction
- White coating
- Anatomically shaped seat and backrest
- Washable synthetic leather, can be disinfected
- 2 rear folding handles
- 2 front telescoping handles
- Folding armrests
- Folding footrest
- 2 stationary and 2 swivel wheels with brakes
- Equipped with 2 pins for locking into floor mount utila model BH 501



	<i>Length Handles Extended</i>	<i>Length Handles Retracted</i>	<i>Width</i>	<i>Height</i>	<i>Seat Height</i>	<i>Weight</i>	<i>Load Limit</i>
<i>Metric</i>	144	75	56	113	48	15	150