



Compact 1, 1-S, 2

The Compact Carry Chair is durable, robust and is manufactured from stainless steel. This proven chair is invaluable when transferring patients through corridors and lifts.

- The one-piece back and seat panel allows the patient to adopt a comfortable position
- Covered in a flame retardant, vinyl coated nylon, it can simply be wiped clean for maximum hygiene
- The design incorporates safety rings which slide over the folding mechanism
- Folds neatly into a compact size for easy stowage
- 7" / 17.8 cm wheels for ease of maneuvering
- The Compact 1-S has the same features as the Compact 1 but features swivel wheels for ease of manoeuvrability
- Developed from the popular Compact Chair, the Compact 2 is manufactured from stainless steel for durability and strength and has internal reinforcements for added rigidity
- For space efficiency the Compact 1 and Compact 2 may be mounted on a wall with optional wall bracket
- The Compact is also available with an optional zipped carry case



Compact 1



Compact 1-S



Compact 2



	Height	Width	Depth	Weight	Load Capacity
Metric	93 cm	44½ cm	42 cm	9½ kg	200 kg

Compact 3



The compact 3 chair has been designed to transport a patient, in the seated position, up and down stairs and also on level ground.

- Innovative tri-wheel system enables smooth operation up and down stairs and reduces requirement for lifting
- The upper carrying handle is fully adjustable allowing easy operation and help to improve operator control
- Folds for compact storage



	<i>Height</i>	<i>Width</i>	<i>Depth</i>	<i>Weight</i>	<i>Load Capacity</i>
Metric	98 cm	41 cm	43 cm	12½ kg	127 kg

Sirocco®

The Sirocco® Easy Transfer System turns the difficult task of patient transportation into an absolute breeze. Lightweight at just 13.5 kg yet robust enough to handle most patient transfer or evacuation situations, Sirocco® can handle up to 159 kg. Operator and patient confidence are enhanced by intelligent, award winning design features.

- Innovative tri-wheel assembly facilitates negotiation of stairs
- Wheel retarding system controls the speed of descent down stairs
- Optional-adjustable, padded head pad ensures patient comfort even if unconscious



	<i>Height</i>	<i>Width</i>	<i>Depth</i>	<i>Weight</i>	<i>Load Capacity</i>
Metric	148 cm	50 cm	96 cm	13½ kg	159 kg

General Evacuation Chair

This easy transfer chair is the ultimate carrying and evacuation chair, combining maximum patient comfort and security with an advanced tri-wheel system, along with fully adjustable foam handle grips at the foot and head ends, to climb and descend stairs with ease and control.

The chair has been specially constructed to handle the most difficult industrial evacuation and rescue situations and can handle patient loads up to 159 kg, while its strong, lightweight construction makes it easy to handle and manoeuvre.

Measures in cm, 148 l x 50 w x 96 d - Weight: 13½ kg.

