


Neil Robertson Stretcher



The Neil Robertson Stretcher is an industry standard rescue stretcher, specifically designed for use in difficult or unconventional rescue circumstances, where a patient needs to be lowered or lifted to safety. Constructed from strong spliced slats and tough cotton, it ensures the patient is secure and stable before and during lifting. Ideal for use in mines, factories or construction sites.



Paraguard™

Paraguard™ has been the official NATO rescue stretcher since WWII, and is used extensively in underground, earthquake and other emergency situations. As it can be used for horizontal and vertical lifting, it is the ideal stretcher for underground, confined space, and rescue from height. Both lightweight and compact, the Paraguard™ is not only easy to carry but easy to use and can be ready in seconds, folding in half making it easy to carry and store.

Other features include:

- Made from high quality aluminium and flame-resistant materials resistant to corrosion and rotting
- Can be carried using four integral side handles that clip out of the way when not in use
- Unique construction allows the centre to unlock and bend for use in awkward spaces
- Impervious to water, grease and oils

