

Saver Model S-2032 Stacking Catastrophe Stretcher




The Saver S-2032 is the ideal product for mass casualty and catastrophic situations. Constructed of durable, lightweight aluminium, the S-2032 conveniently and quickly deploys onto the rescue scene to accommodate multiple victims. Once finished, it conveniently stacks onto itself to provide the most compact storage available. The heavy duty vinyl-coated nylon cover is durable, fire retardant, and stain resistant to blood and bodily fluids.



S2032 Stackable Catastrophe Stretcher

Length	187 cm
Width	48 cm
Height	18 cm
Unit Weight	3 kg
Loading Capacity	159 kg



Saver Model S-2033 XWF Stretcher

The Saver S-2033 XWF Stretcher is an extra wide and foldable stretcher with two I.V. pole holders allowing it to function as a field bed in emergency situations.

- Extra wide I.V. pole holder on each side
- Foldable (secured by locking device) for easier storage and transport

0153216 Model S-2033, Xtra Wide Foldable Stretcher

Length	195 cm
Width	80 cm
Height	45 cm
Unit Weight	10 kg
Loading Capacity	159 kg





Saver Model S-2108 Emergency Pole Stretcher Series

The Saver Folding Pole Stretcher series is simple in design, easy-to-use, and durable enough to get the job done. It features heavy-duty aluminum poles for easy handling. The rugged steel spreader bars and aluminium feet are standard. A heavy-duty vinyl-coated nylon cover is durable, fire retardant, and stain resistant while the no-slip handgrips help technicians keep a safe grip on the heaviest of loads. The A-style collapses into itself crosswise, while the AF version collapses both crosswise and in half lengthwise for the most compact storage available. Also available in triple-folding version or with retractable handles.



0153166 Model S-2108-A, Pole Stretcher Series

Length	206 cm
Folded	229 cm
Width: Open	55 cm
Width : Folded	11 cm
Height	14 cm
Unit Weight	5 kg
Loading Capacity	159 kg

0153168 Model S-2108-AF, Pole Stretcher Series

Length	206/229 cm
Folded	103/114 cm
Width: Open	55 cm
Width : Folded	11 cm
Height	20 cm
Unit Weight	6 kg
Loading Capacity	159 kg



Saver Model S-2107 Combination Stretcher Chair

Designed to serve as three patient handling devices in one: a wheeled chair, a stair chair, and a flat stretcher. Sturdy, light-weight aluminium frame provides long-term durability, while the vinyl-coated nylon cover wears well and is easy to clean and maintain. Folds in half for compact storage. Features a highback chair design that supports the patient's head so personnel can maintain proper posture when moving the patient. 10 cm wheels enable easy movement over most surfaces. Available in burgundy and orange.



0153221 Model S-2107, Combination Stretcher Chair

Length	185 cm
Folded	97 cm
Width: Open	49 cm
Width : Folded	55 cm
Height	20 cm
Unit Weight	5 kg
Loading Capacity	159 kg

Includes two Saver 210 one-piece Black restraints with plastic buckles, and one Saver 125 one-piece Black restraint with plastic buckles.





Saver Model S-212 Folding Emergency Stretcher

The Saver S-212 is a simple, flat emergency stretcher constructed of durable, lightweight aluminium. The stretcher folds easily for quick and compact storage, and comes with two patient restraints. The heavy-duty vinyl-coated nylon cover is durable, fire retardant, and stain resistant to blood and bodily fluids.

0153171 Model S-212, Folding Emergency Stretcher

Length	187 cm
Folded	93 cm
Width:	48 cm
Height : Open	3 cm
Height : Folded	6 cm
Unit Weight	4 kg
Loading Capacity	159 kg

Includes two Saver150 Black two-piece restraints with plastic buckles.



Saver Model S-211 Folding Emergency Stretcher

The Saver S-211 is a simple, flat emergency stretcher constructed of durable, lightweight aluminium. The stretcher folds easily for quick and compact storage, and comes with two patient restraints. The heavy-duty vinyl-coated nylon cover is durable, fire retardant, and stain resistant to blood and bodily fluids. The Model S-211 features two 10 cm ground wheels and two stationary feet to make patient transfer easy.

0153170 Model S-211, Folding Emergency Stretcher

Length	187 cm
Folded	93 cm
Width:	48 cm
Height : Open	11 cm
Height : Folded	16 cm
Unit Weight	6 kg
Loading Capacity	159 kg

Includes two Saver150 Black two-piece restraints with plastic buckles.





Saver Model S-29 Emergency Stretcher with wheels & backrest

- Designed for vehicles with limited head room
- Head-end frame section telescopes to shorten stretcher
- Backrest is adjustable to eight angles (0° to 75°)
- 10 cm wheels at head-end with stationary posts at foot-end
- Folds in half for compact storage

0153220 Model S-29, Emergency Stretcher with wheels and backrest

Length	191½ cm
Folded	93 cm
Width:	49 cm
Height : Open	15½ cm
Height : Folded	20½ cm
Unit Weight	8 kg
Loading Capacity	159 kg



Includes two Saver150 Black two-piece restraints with plastic buckles.



Saver Model S-2409 Lift-Off Stretcher

The Saver S-2409 Lift-off Stretcher adjusts to a variety of positions and is multi-functional, designed to accommodate a wide range of patient needs and emergency situations. Stretcher may be reversed to enable patient to be loaded head-first or feet-first.

- Adjustable backrest mechanism
- Swing-down side rails enable convenient patient transfer from cot to bed
- Ask for undercarriage compatibility



0153222 Model S-2409, Lift-Off Stretcher

Length	187 cm
Width	58 cm
Height	20-90 cm
Unit Weight	21 kg
Loading Capacity	159 kg



Includes two Saver150 Black two-piece restraints with plastic buckles and orange segmented mattress.

