



Saver Model S-219 Cot

The Saver S-219 is easy to use and built tough for years of reliable service. The adjustable backrest can be easily locked in any of eight positions, up to an angle of 76 degrees, providing increased comfort to the recumbent patient. Swivel wheels allow smooth and swift movement in any direction. Two of the wheels include foot brakes for added safety.

0054159 Model S-219, Cot

Length	238½ cm
Width:	51½ cm
Height : Open	110 cm
Height : Folded	41 cm
Unit Weight	26 kg
Loading Capacity	159 kg

Available with two Saver 150 black two-piece restraints with metal buckles and mattress.



Saver Model S-226 Cot

The Saver S-226 Roll-in Cot has two levels of height adjustment: fully elevated and ground level. Adjustable backrest angle from 0° to 75°. Features an automatic leg locking system that ensures cot legs are locked in place once they are fully unfolded. Optional ambulance fastener with 10 G cot lock.

0054161 Model S-226, Cot

Length	187 cm
Folded	93 cm
Width:	48 cm
Height : Open	11 cm
Height : Folded	16 cm
Unit Weight	6 kg
Loading Capacity	159 kg

