



Saver Model S-210 Scout Carry Chair

- Aluminium construction, folds for easy storage, and it is ideal for small, confined areas.
- Two rear wheels (7.6 cm) for easy rolling
- Includes 2 one-piece restraints
- Folds for compact storage
- Sewn orange vinyl cover

0721326 Model S-210, Scout Carry Chair

Length	57½ cm
Folded	17cm
Width:	45cm
Height : Open	96cm
Height : Folded	10¾ cm
Unit Weight	3 kg
Loading Capacity	159 kg

Includes two Saver210 one - piece black restraints with plastic buckles.



Saver Model S-214 Patient Carry Chair

- Stainless steel construction
- Folds for easy storage
- Dual purpose handles: patient comfort and operator control
- Four Wheels: two stationary (rear) and two Swivel (front)
- Sewn orange vinyl cover



0721327 Model S-214, Patient Carry Chair

Length (handles up)	74½ cm
Length (handles down)	82.2 cm
Length (folded)	50.8 cm
Width (handles up)	51.4 cm
Width (handles down)	54 cm
Width (folded)	78.7 cm
Height (handles up)	94 cm
Height (handles down)	89.8 cm
Height (folded)	19.7 cm
Unit Weight	10.8 cm
Loading Capacity	159 kg





Saver Model S-240 Folding Patient Transport Chair

Easy to maneuver by two personnel, the Saver S-240 Patient Transport Chair gets you in and out of cramped quarters. Two 12.7 cm stationary rear wheels enable personnel to move the patient over most surfaces with minimal effort. Telescoping, two-position front handles allow foot-end personnel to maneuver the patient without interfering with the patient's feet. Vinyl-coated nylon cover is simple to clean, resists stains, blood and body fluids. Includes three patient restraints. Optional vehicle fastener.



0721324 Model S-240, Folding Chair

Length	121 cm
Folded	17 cm
Width:	52 cm
Height : Open	91½ cm
Height : Folded	91½ cm
Unit Weight	9 kg
Loading Capacity	159 kg

Includes two Saver180 one-piece black restraints with plastic buckles and one Saver 50 mm one-piece black restraint with plastic buckles.



Saver Model S-242 Folding Patient Transport Chair

Easy to maneuver by two personnel, the Saver S-242 Patient Transport Chair gets you in and out of cramped quarters. Two 12.7 cm stationary rear wheels and two 7.6 cm, front mounted swivel wheels enable personnel to move the patient over most surfaces with minimal effort. Telescoping, two-position front handles allow foot-end personnel to maneuver the patient without interfering with the patient's feet. Vinyl-coated nylon cover is simple to clean, resists stains, blood and body fluids. Includes three patient restraints. Vehicle fastener available.



0721323 Model S-242, Folding Patient Transport Chair

Length	121 cm
Folded	17 cm
Width:	52 cm
Height : Open	91½ cm
Height : Folded	91½ cm
Unit Weight	10 kg
Loading Capacity	159 kg

Includes two Saver180 one-piece black restraints with plastic buckles and one Saver 50 mm one-piece black restraint with plastic buckles.

