



Saver Model S-265 Scoop Stretcher

Designed so EMS personnel can uncouple either end or both ends of the stretcher and gently scoop up the patient using scissor-type closing motion. Ideal for hip injuries and accident victims found on the street. Supports a patient in the position found, reducing risk of further injury. Features include a narrow foot-end frame for handling in confined areas, telescoping tubes allow length adjustments to fit patients of various heights, open-center design permits patient to be X-rayed.

- Twin Safety Lock for smooth locking and unlocking
- Blue painted surface
- Includes three patient restraints
- Optional vehicle fastener available



0153219 Model S-265, Scoop Stretcher

Length	166-202 cm
Folded	120 cm
Width:	42.2 cm
Height : Open	6.7 cm
Height : Folded	10 cm
Unit Weight	9¾ kg
Loading Capacity	159 kg

Includes three Saver 210 one - piece black restraints with plastic buckles.

