

Saver KED350



- Improved means of immobilizing and extricating patients from auto accidents or confined spaces
- The wrap-around design provides horizontal flexibility for application, and vertical rigidity for support
- Immobilizes torso, head and neck minimizing risks of further injury
- Color coded straps for visibility
- X-rays or advanced life support procedures can be performed with the KED350 in place
- Heavy-duty, vinyl-coated nylon is durable and easy to clean
- Case is large enough to carry extrication collars
- Built-in handles enable rescuers to get a firm grip on the patient and KED350 as one unit
- Can be inverted as hip splint; adapts for children and pregnant women
- Adjustable fold-back sides permit easy access to patient's chest
- Can be used for multiple patient sizes, including pediatric patients, children, and pregnant women
- Stores rolled for immediate use
- Available in blue only
- Sold in singles or pack of 8

0314022 KED350, one each

000KED350 KED350, pack of 8



Load Limit: 159 kg (350 lb)

Shipping Dimensions (pack of 8): 86x53x38 cm (34"x21"x15")

Shipping Weight (pack of 8): 26 kg (57 lb)



Visit SAVER website www.saver-ems.com for more information.

To order SAVER products, contact your Ferno Dealer,

or E-mail: internationalsales@saver-ems.com

email: internationalsales@saver-ems.com